

Androgen-dependent TFPI-regulating protein

Catalog No: #AP80169



Package Size: #AP80169-1 50ug #AP80169-2 100ug #AP80169-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Androgen-dependent TFPI-regulating protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Liver regeneration-related protein LRRG140,Aa2-020
Accession No.	Uniprot ID:Q5M828
Uniprot	Q5M828
GeneID	361228;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adtrp

Note: This product is for in vitro research use only