

TNF α Monoclonal antibody(Q36)

Catalog No: #44074

Package Size: #44074-1 50ul #44074-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	TNF α Monoclonal antibody(Q36)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q36
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	The antibody detects endogenous TNF α protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	TNF α
Other Names	TNF α
Accession No.	Swiss-Prot:P01375Gene ID:7124
Uniprot	P01375
GeneID	7124;
Calculated MW	17,25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:1000-1:3000

Immunohistochemistry: 1:50-1:200

Images

Western blot analysis of Recombinant Human TNF α Protein with TNF α Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using TNF a Mouse mAb diluted at 1:50.

Note: This product is for in vitro research use only