

Disintegrin and metalloproteinase domain-containing protein 28



Catalog No: #AP80172

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80172-1 50ug #AP80172-2 100ug #AP80172-3 1mg

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 28
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM 28, Epididymal metalloproteinase-like, disintegrin-like, and cysteine-rich protein II, eMDC II, Metalloproteinase-like, disintegrin-like, and cysteine-rich protein L, MDC-L, ADAM23, MDCL
Accession No.	Uniprot ID: Q9UKQ2
Uniprot	Q9UKQ2
GeneID	10863;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADAM28

Note: This product is for in vitro research use only