

ADP/ATP translocase 2

Catalog No: #AP80173



Package Size: #AP80173-1 50ug #AP80173-2 100ug #AP80173-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADP/ATP translocase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ADP,ATP carrier protein 2,Adenine nucleotide translocator 2,ANT 2,Solute carrier family 25 member 5,ANT2
Accession No.	Uniprot ID:Q8SQH5
Uniprot	Q8SQH5
GeneID	282479;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC25A5

Note: This product is for in vitro research use only