

Adipocyte enhancer-binding protein 1

Catalog No: #AP80174



Package Size: #AP80174-1 50ug #AP80174-2 100ug #AP80174-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adipocyte enhancer-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AE-binding protein 1,Aortic carboxypeptidase-like protein,ACL P
Accession No.	Uniprot ID:Q8IUX7
Uniprot	Q8IUX7
GeneID	165;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AE BP1

Note: This product is for in vitro research use only