

Aminoacyl tRNA synthase complex-interacting multifunctional protein 1

Catalog No: #AP80175

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80175-1 50ug #AP80175-2 100ug #AP80175-3 1mg

Description

| | |
|--------------------|---|
| Product Name | Aminoacyl tRNA synthase complex-interacting multifunctional protein 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Multisynthase complex auxiliary component p43, Emap2, Scye1 |
| Accession No. | Uniprot ID:P31230 |
| Uniprot | P31230 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Aimp1

Note: This product is for in vitro research use only