

Acyl-CoA desaturase

Catalog No: #AP80177



Package Size: #AP80177-1 50ug #AP80177-2 100ug #AP80177-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Acyl-CoA desaturase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Delta(9)-desaturase,Delta-9 desaturase,Fatty acid desaturase,Stearoyl-CoA desaturase |
| Accession No. | Uniprot ID:Q9TT94 |
| Uniprot | Q9TT94 |
| GeneID | 280924; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SCD

Note: This product is for in vitro research use only