

Aspartoacylase

Catalog No: #AP80178



Package Size: #AP80178-1 50ug #AP80178-2 100ug #AP80178-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Aspartoacylase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Aminoacylase-2,ACY-2,ACY2 |
| Accession No. | Uniprot ID:P46446 |
| Uniprot | P46446 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ASPA

Note: This product is for in vitro research use only