

P53(Mono Methyl Lys370) Monoclonal antibody(Q59)

Catalog No: #HW190

Package Size: #HW190-1 50ul #HW190-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	P53(Mono Methyl Lys370) Monoclonal antibody(Q59)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q59
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	The antibody detects endogenous P53 (Mono Methyl Lys370) protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	P53
Other Names	P53(Mono Methyl Lys370)
Accession No.	Swiss-Prot:P04637Gene ID:7157
Uniprot	P04637
GeneID	7157;
Calculated MW	53kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using P53(Mono Methyl Lys370) Mouse mAb diluted at 1:500.

Note: This product is for in vitro research use only