

7-methylguanosine phosphate-specific 5'-nucleotidase

Catalog No: #AP80180



Package Size: #AP80180-1 50ug #AP80180-2 100ug #AP80180-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	7-methylguanosine phosphate-specific 5'-nucleotidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	7-methylguanosine nucleotidase,Cytosolic 5'-nucleotidase 3B,Cytosolic 5'-nucleotidase III-like protein,cN-III-like protein,N(7)-methylguanylate 5'-phosphatase,NT5C3L
Accession No.	Uniprot ID:Q969T7
Uniprot	Q969T7
GeneID	115024;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NT5C3B

Note: This product is for in vitro research use only