

## 5-hydroxytryptamine receptor 2A

Catalog No: #AP80182



Package Size: #AP80182-1 50ug #AP80182-2 100ug #AP80182-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | 5-hydroxytryptamine receptor 2A   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Dog   |
| Other Names        | 5-HT-2,5-HT-2A,5-hydroxytryptamine receptor,Serotonin receptor 2A,5-HTR2A |
| Accession No.      | Uniprot ID:O46635   |
| Uniprot            | O46635  |
| GeneID             | 403882;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                    |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HTR2A

Note: This product is for in vitro research use only