

# Amyloid beta A4 protein

Catalog No: #AP80183



Package Size: #AP80183-1 50ug #AP80183-2 100ug #AP80183-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Amyloid beta A4 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	ABPP,APP,Alzheimer disease amyloid A4 protein homolog
Accession No.	Uniprot ID:Q28280
Uniprot	Q28280
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:APP

Note: This product is for in vitro research use only