

Alpha-actinin-4

Catalog No: #AP80188



Package Size: #AP80188-1 50ug #AP80188-2 100ug #AP80188-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-actinin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Non-muscle alpha-actinin 4
Accession No.	Uniprot ID:A5D7D1
Uniprot	A5D7D1
GeneID	522269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ACTN4

Note: This product is for in vitro research use only