

Beta-1 adrenergic receptor

Catalog No: #AP80190



Package Size: #AP80190-1 50ug #AP80190-2 100ug #AP80190-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-1 adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-1 adrenoreceptor, Beta-1 adrenoceptor, Adrb1r
Accession No.	Uniprot ID:P18090
Uniprot	P18090
GeneID	24925;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adrb1

Note: This product is for in vitro research use only