

Neuronal acetylcholine receptor subunit alpha-4

Catalog No: #AP80203



Package Size: #AP80203-1 50ug #AP80203-2 100ug #AP80203-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neuronal acetylcholine receptor subunit alpha-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Acra4
Accession No.	Uniprot ID:P09483
Uniprot	P09483
GeneID	25590;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Chrna4

Note: This product is for in vitro research use only