

Adenosine deaminase-like protein

Catalog No: #AP80212

Package Size: #AP80212-1 50ug #AP80212-2 100ug #AP80212-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine deaminase-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q80SY6
Uniprot	Q80SY6
GeneID	75894;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adal

Note: This product is for in vitro research use only