

ADP/ATP translocase 1

Catalog No: #AP80215



Package Size: #AP80215-1 50ug #AP80215-2 100ug #AP80215-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADP/ATP translocase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADP,ATP carrier protein 1,ADP,ATP carrier protein, heart/skeletal muscle isoform T1,Adenine nucleotide translocator 1,ANT 1,Solute carrier family 25 member 4,mANC1,Anc1,Ant1
Accession No.	Uniprot ID:P48962
Uniprot	P48962
GeneID	11739;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc25a4

Note: This product is for in vitro research use only