

Type-1B angiotensin II receptor

Catalog No: #AP80216



Package Size: #AP80216-1 50ug #AP80216-2 100ug #AP80216-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Type-1B angiotensin II receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AT3,Angiotensin II type-1B receptor,AT1B,Agtr1,At1b
Accession No.	Uniprot ID:P29089
Uniprot	P29089
GeneID	81638;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Agtr1b

Note: This product is for in vitro research use only