

Amiloride-sensitive amine oxidase [copper-containing]

Catalog No: #AP80217



Package Size: #AP80217-1 50ug #AP80217-2 100ug #AP80217-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amiloride-sensitive amine oxidase [copper-containing]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DAO,Diamine oxidase,Amiloride-binding protein,ABP,Histaminase
Accession No.	Uniprot ID:P36633
Uniprot	P36633
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abp1

Note: This product is for in vitro research use only