

Neuronal acetylcholine receptor subunit alpha-2

Catalog No: #AP80218



Package Size: #AP80218-1 50ug #AP80218-2 100ug #AP80218-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuronal acetylcholine receptor subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Acra2
Accession No.	Uniprot ID:P12389
Uniprot	P12389
GeneID	170945;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Chrna2

Note: This product is for in vitro research use only