

Disintegrin and metalloproteinase domain-containing protein 17

Catalog No: #AP80220

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80220-1 50ug #AP80220-2 100ug #AP80220-3 1mg

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM 17, TNF-alpha convertase, TNF-alpha-converting enzyme, CD156b, Tace
Accession No.	Uniprot ID: Q9Z0F8
Uniprot	Q9Z0F8
GeneID	11491;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adam17

Note: This product is for in vitro research use only