

Angiotensin-converting enzyme

Catalog No: #AP80225



Package Size: #AP80225-1 50ug #AP80225-2 100ug #AP80225-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ACE,Dipeptidyl carboxypeptidase I,Kininase II,CD143,Dcp1
Accession No.	Uniprot ID:P47820
Uniprot	P47820
GeneID	24310;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ace

Note: This product is for in vitro research use only