

5'-AMP-activated protein kinase catalytic subunit alpha-1

Catalog No: #AP80233

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80233-1 50ug #AP80233-2 100ug #AP80233-3 1mg

Description

Product Name	5'-AMP-activated protein kinase catalytic subunit alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	AMPK subunit alpha-1,AMPK 63 kDa subunit,Acetyl-CoA carboxylase kinase,ACACA kinase,Hydroxymethylglutaryl-CoA reductase kinase,HMGCR kinase,Tau-protein kinase PRKAA1,AMPK1
Accession No.	Uniprot ID:Q09136
Uniprot	Q09136
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKAA1

Note: This product is for in vitro research use only