

Adenylate cyclase type 1

Catalog No: #AP80234



Package Size: #AP80234-1 50ug #AP80234-2 100ug #AP80234-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenylate cyclase type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ATP pyrophosphate-lyase 1, Adenylate cyclase type I, Adenylyl cyclase 1, AC1, Ca(2+)/calmodulin-activated adenylyl cyclase
Accession No.	Uniprot ID:P19754
Uniprot	P19754
GeneID	281601;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADCY1

Note: This product is for in vitro research use only