

Aconitate hydratase, mitochondrial

Catalog No: #AP80235



Package Size: #AP80235-1 50ug #AP80235-2 100ug #AP80235-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aconitate hydratase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Aconitase,Citrate hydro-lyase
Accession No.	Uniprot ID:Q99K10
Uniprot	Q99K10
GeneID	11429;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aco2

Note: This product is for in vitro research use only