

Alpha-galactosidase A

Catalog No: #AP80239



Package Size: #AP80239-1 50ug #AP80239-2 100ug #AP80239-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-galactosidase A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-D-galactosidase A,Alpha-D-galactoside galactohydrolase,Melibiase,Ags
Accession No.	Uniprot ID:P51569
Uniprot	P51569
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gla

Note: This product is for in vitro research use only