

Medium-chain specific acyl-CoA dehydrogenase, mitochondrial

Catalog No: #AP80240

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80240-1 50ug #AP80240-2 100ug #AP80240-3 1mg

Description

Product Name	Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MCAD
Accession No.	Uniprot ID:Q3SZB4
Uniprot	Q3SZB4
GeneID	505968;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ACADM

Note: This product is for in vitro research use only