

Alpha-2A adrenergic receptor

Catalog No: #AP80243



Package Size: #AP80243-1 50ug #AP80243-2 100ug #AP80243-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-2A adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Alpha-2A adrenoreceptor,Alpha-2A adrenoceptor,Alpha-2D adrenergic receptor
Accession No.	Uniprot ID:Q28838
Uniprot	Q28838
GeneID	282135;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADRA2A

Note: This product is for in vitro research use only