

Single-stranded DNA cytosine deaminase

Catalog No: #AP80244



Package Size: #AP80244-1 50ug #AP80244-2 100ug #AP80244-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Single-stranded DNA cytosine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activation-induced cytidine deaminase,AID,Cytidine aminohydrolase,AID
Accession No.	Uniprot ID:Q9GZX7
Uniprot	Q9GZX7
GeneID	57379;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AICDA

Note: This product is for in vitro research use only