

5-hydroxytryptamine receptor 2A

Catalog No: #AP80247



Package Size: #AP80247-1 50ug #AP80247-2 100ug #AP80247-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5-hydroxytryptamine receptor 2A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	5-HT-2,5-HT-2A,5-hydroxytryptamine receptor,Serotonin receptor 2A,5-HTR2A
Accession No.	Uniprot ID:O46635
Uniprot	O46635
GeneID	403882;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HTR2A

Note: This product is for in vitro research use only