

Cytosolic 5'-nucleotidase III-like protein

Catalog No: #AP80251



Package Size: #AP80251-1 50ug #AP80251-2 100ug #AP80251-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosolic 5'-nucleotidase III-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	cN-III-like protein,RCJMB04_11c18
Accession No.	Uniprot ID:Q5ZKF6
Uniprot	Q5ZKF6
GeneID	420034;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NT5C3L

Note: This product is for in vitro research use only