

Amiloride-sensitive amine oxidase [copper-containing]

Catalog No: #AP80253



Package Size: #AP80253-1 50ug #AP80253-2 100ug #AP80253-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amiloride-sensitive amine oxidase [copper-containing]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	DAO, Diamine oxidase, Amiloride-binding protein, ABP, Histaminase, PK-DAO
Accession No.	Uniprot ID:Q9TRC7
Uniprot	Q9TRC7
GeneID	100517436;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ABP1

Note: This product is for in vitro research use only