

Angiotensin-converting enzyme 1

Catalog No: #AP80257

Package Size: #AP80257-1 50ug #AP80257-2 100ug #AP80257-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Angiotensin-converting enzyme 1
Accession No.	Uniprot ID:Q1RMR1
Uniprot	Q1RMR1
GeneID	509971;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ANGPTL1

Note: This product is for in vitro research use only