

AP-2 complex subunit beta

Catalog No: #AP80258



Package Size: #AP80258-1 50ug #AP80258-2 100ug #AP80258-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AP-2 complex subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AP105B,Adaptor protein complex AP-2 subunit beta,Adaptor-related protein complex 2 subunit beta,Beta-2-adaptin,Beta-adaptin,Clathrin assembly protein complex 2 beta large chain,Plasma membrane adaptor HA2/AP2 adaptin beta subunit,ADTB2,CLAPB1
Accession No.	Uniprot ID:P63010
Uniprot	P63010
GeneID	163;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AP2B1

Note: This product is for in vitro research use only