

ATP-binding cassette sub-family G member 2

Catalog No: #AP80259

Package Size: #AP80259-1 50ug #AP80259-2 100ug #AP80259-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-binding cassette sub-family G member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Brain multidrug resistance protein,Urate exporter,CD338,BMDP
Accession No.	Uniprot ID:Q8MIB3
Uniprot	Q8MIB3
GeneID	397073;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ABCG2

Note: This product is for in vitro research use only