

Disintegrin and metalloproteinase domain-containing protein 9

Catalog No: #AP80261

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80261-1 50ug #AP80261-2 100ug #AP80261-3 1mg

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM 9, Cellular disintegrin-related protein, Meltrin-gamma, Metalloprotease/disintegrin/cysteine-rich protein 9, Myeloma cell metalloproteinase, KIAA0021, MCMP, MDC9, MLTNG
Accession No.	Uniprot ID:Q13443
Uniprot	Q13443
GeneID	8754;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADAM9

Note: This product is for in vitro research use only