

Allograft inflammatory factor 1

Catalog No: #AP80262



Package Size: #AP80262-1 50ug #AP80262-2 100ug #AP80262-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Allograft inflammatory factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AIF-1,Ionized calcium-binding adapter molecule 1,lba1
Accession No.	Uniprot ID:O70200
Uniprot	O70200
GeneID	11629;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aif1

Note: This product is for in vitro research use only