

A-kinase anchor protein 1, mitochondrial

Catalog No: #AP80263



Package Size: #AP80263-1 50ug #AP80263-2 100ug #AP80263-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	A-kinase anchor protein 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	A-kinase anchor protein 121 kDa,AKAP 121,Dual specificity A-kinase-anchoring protein 1,D-AKAP-1,Protein kinase A-anchoring protein 1,PRKA1,Spermatid A-kinase anchor protein 84,S-AKAP84,Akap
Accession No.	Uniprot ID:O88884
Uniprot	O88884
GeneID	114124;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Akap1

Note: This product is for in vitro research use only