

ATP-citrate synthase

Catalog No: #AP80264



Package Size: #AP80264-1 50ug #AP80264-2 100ug #AP80264-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ATP-citrate synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ATP-citrate (pro-S-)-lyase,Citrate cleavage enzyme
Accession No.	Uniprot ID:Q32PF2
Uniprot	Q32PF2
GeneID	511135;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ACLY

Note: This product is for in vitro research use only