

## ADP/ATP translocase 1

Catalog No: #AP80266



Package Size: #AP80266-1 50ug #AP80266-2 100ug #AP80266-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ADP/ATP translocase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ADP,ATP carrier protein 1,ADP,ATP carrier protein, heart isoform T1,Adenine nucleotide translocator 1,ANT 1,Solute carrier family 25 member 4,ANT1
Accession No.	Uniprot ID:P02722
Uniprot	P02722
GeneID	282478;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC25A4

Note: This product is for in vitro research use only