

Amiloride-sensitive amine oxidase [copper-containing]

Catalog No: #AP80268



Package Size: #AP80268-1 50ug #AP80268-2 100ug #AP80268-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amiloride-sensitive amine oxidase [copper-containing]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DAO,Diamine oxidase,Amiloride-binding protein,ABP,Histaminase,Kidney amine oxidase,KAO,AOC1,DAO1
Accession No.	Uniprot ID:P19801
Uniprot	P19801
GeneID	26;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ABP1

Note: This product is for in vitro research use only