

Amiloride-sensitive amine oxidase [copper-containing]

Catalog No: #AP80269



Package Size: #AP80269-1 50ug #AP80269-2 100ug #AP80269-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------------------------------------------|
| Product Name | Amiloride-sensitive amine oxidase [copper-containing] |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | DAO,Diamine oxidase,Amiloride-binding protein,ABP,Histaminase |
| Accession No. | Uniprot ID:Q8JZQ5 |
| Uniprot | Q8JZQ5 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abp1

Note: This product is for in vitro research use only