

Agouti-related protein

Catalog No: #AP80270



Package Size: #AP80270-1 50ug #AP80270-2 100ug #AP80270-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Agouti-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AGRT,ART
Accession No.	Uniprot ID:P56413
Uniprot	P56413
GeneID	280987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AGRP

Note: This product is for in vitro research use only