

Adenosine deaminase domain-containing protein 1

Catalog No: #AP80273



Package Size: #AP80273-1 50ug #AP80273-2 100ug #AP80273-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine deaminase domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Testis nuclear RNA-binding protein, TENR
Accession No.	Uniprot ID:Q96M93
Uniprot	Q96M93
GeneID	132612;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADAD1

Note: This product is for in vitro research use only