

Adenylate cyclase type 8

Catalog No: #AP80274



Package Size: #AP80274-1 50ug #AP80274-2 100ug #AP80274-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenylate cyclase type 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATP pyrophosphate-lyase 8, Adenylate cyclase type VIII, Adenylyl cyclase 8, AC8, Ca(2+)/calmodulin-activated adenylyl cyclase
Accession No.	Uniprot ID:P40145
Uniprot	P40145
GeneID	114;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: ADCY8

Note: This product is for in vitro research use only