

Alpha-actinin-3

Catalog No: #AP80279



Package Size: #AP80279-1 50ug #AP80279-2 100ug #AP80279-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-actinin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-actinin skeletal muscle isoform 3,F-actin cross-linking protein
Accession No.	Uniprot ID:O88990
Uniprot	O88990
GeneID	11474;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Actn3

Note: This product is for in vitro research use only