

Natriuretic peptides A

Catalog No: #AP80281



Package Size: #AP80281-1 50ug #AP80281-2 100ug #AP80281-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Natriuretic peptides A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Prepronatriodilatin
Accession No.	Uniprot ID:P07501
Uniprot	P07501
GeneID	281355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NPPA

Note: This product is for in vitro research use only