

Membrane primary amine oxidase

Catalog No: #AP80285



Package Size: #AP80285-1 50ug #AP80285-2 100ug #AP80285-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Membrane primary amine oxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Copper amine oxidase,Semicarbazide-sensitive amine oxidase,SSAO,Vascular adhesion protein 1,VAP-1,Vap1
Accession No.	Uniprot ID:O70423
Uniprot	O70423
GeneID	11754;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aoc3

Note: This product is for in vitro research use only