

Aldo-keto reductase family 1 member C1

Catalog No: #AP80286



Package Size: #AP80286-1 50ug #AP80286-2 100ug #AP80286-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Aldo-keto reductase family 1 member C1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | 20-alpha-hydroxysteroid dehydrogenase,20-alpha-HSD,Chlordecone reductase homolog HAKRC,Dihydrodiol dehydrogenase 1/2,DD1/DD2,High-affinity hepatic bile acid-binding protein,HBAB,Indanol dehydrogenase,Trans-1,2-dihydrobenzene-1,2-diol dehydrogenase,DDH,DDH1 |
| Accession No. | Uniprot ID:Q04828 |
| Uniprot | Q04828 |
| GeneID | 1645; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AKR1C1

Note: This product is for in vitro research use only