

Alpha-1A adrenergic receptor

Catalog No: #AP80291



Package Size: #AP80291-1 50ug #AP80291-2 100ug #AP80291-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-1A adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alpha-1A adrenoreceptor,Alpha-1A adrenoceptor,Alpha-1C adrenergic receptor,Alpha-adrenergic receptor 1c,ADRA1C
Accession No.	Uniprot ID:P35348
Uniprot	P35348
GeneID	148;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADRA1A

Note: This product is for in vitro research use only